Package 'ptw'

January 19, 2022

Type Package

Title Parametric Time Warping

Version 1.9-16

Description Parametric Time Warping aligns patterns, i.e. it aims to put corresponding features at the same locations. The algorithm searches for an optimal polynomial describing the warping. It is possible to align one sample to a reference, several samples to the same reference, or several samples to several references. One can choose between calculating individual warpings, or one global warping for a set of samples and one reference. Two optimization criteria are implemented: RMS (Root Mean Square error) and WCC (Weighted Cross Correlation). Both warping of peak profiles and of peak lists are supported. A vignette for the latter is contained in the inst/doc directory of the source package - the vignette source can be found on the package github site.

URL https://github.com/rwehrens/ptw

License GPL (>= 2)

Imports RcppDE, graphics, grDevices, stats

NeedsCompilation yes

Author Jan Gerretzen [ctb], Paul Eilers [aut], Hans Wouters [ctb], Tom Bloemberg [aut], Ron Wehrens [aut, cre]

Maintainer Ron Wehrens <ron.wehrens@gmail.com> Repository CRAN Date/Publication 2022-01-19 13:42:42 UTC

R topics documented:

asysm	3
baseline.corr	4
bestref	5
calc.multicoef	6
calc.zerocoef	7
coda	8
difsm	9
gaschrom	10
lcms	11
mzchannel2pktab	12
padzeros	13
plot.ptw	
predict.ptw	16
ptw	18
ptwgrid	23
RMS	24
select.traces	26
warp.time	27
wcc	28
whit1	29
whit2	30
	31

Index

ptw-package

Parametric Time Warping

Description

Parametric Time Warping aligns patterns, i.e. it aims to put corresponding features at the same locations. The algorithm searches for an optimal polynomial describing the warping. It is possible to align one sample to a reference, several samples to the same reference, or several samples to several references. One can choose between calculating individual warpings, or one global warping for a set of samples and one reference. Two optimization criteria are implemented: RMS (Root Mean Square error) and WCC (Weighted Cross Correlation). Both warping of peak profiles and of peak lists are supported. A vignette for the latter is contained in the inst/doc directory of the source package - the vignette source can be found on the package github site.

Details

The DESCRIPTION file:

Package:	ptw
Type:	Package
Title:	Parametric Time Warping
Version:	1.9-16
Authors@R:	c(person("Jan", "Gerretzen", role = "ctb"), person("Paul", "Eilers", role = "aut"), person("Hans", "Wouters", ro
Description:	Parametric Time Warping aligns patterns, i.e. it aims to put corresponding features at the same locations. The a

asysm

URL:	https://github.com/rwehrens/ptw
License:	GPL (>= 2)
Imports:	RcppDE, graphics, grDevices, stats
Author:	Jan Gerretzen [ctb], Paul Eilers [aut], Hans Wouters [ctb], Tom Bloemberg [aut], Ron Wehrens [aut, cre]
Maintainer:	Ron Wehrens <ron.wehrens@gmail.com></ron.wehrens@gmail.com>

Author(s)

NA

Maintainer: NA

References

@ArticleBloemberg2010, title = Improved Parametric Time Warping for Proteomics, author = Tom G. Bloemberg and Jan Gerretzen and Hans J. P. Wouters and Jolein Gloerich and Maurice van Dael and Hans J. C. T. Wessels and Lambert P. van den Heuvel and Paul H. C. Eilers and Lutgarde M. C. Buydens and Ron Wehrens, year = 2010, journal = Chemometrics and Intelligent Laboratory Systems, volume = 104, number = 1, pages = 65-74,

ArticleWehrens2015, title = Fast parametric warping of peak lists, author = Ron Wehrens and Tom G. Bloemberg and Paul H. C. Eilers, year = 2015, volume = 31, pages = 3063–3065, journal = Bioinformatics, doi = 10.1093/bioinformatics/btv299,

asysm

Trend estimation with asymmetric least squares

Description

Estimates a trend based on asymmetric least squares. In this case used to estimate the baseline of a given spectrum.

Usage

```
asysm(y, lambda = 1e+07, p = 0.001, eps = 1e-8, maxit = 25)
```

Arguments

У	data: either a vector or a data matrix containing spectra as rows
lambda	smoothing parameter (generally 1e5 - 1e8)
р	asymmetry parameter
eps	numerical precision for convergence
maxit	max number of iterations. If no convergence is reached, a warning is issued.

Details

Asymmetric least squares (not to be confused with alternating least squares) assigns different weights to the data points that are above and below an iteratively estimated trendline, respectively. In this case, the asymmetry parameter p ($0 \le p \le 1$) is the weight for points above the trendline, whereas 1-p is the weight for points below it. Naturally, p should be small for estimating baselines. The parameter lambda controls the amount of smoothing: the larger it is, the smoother the trendline will be.

Value

An estimated baseline

Author(s)

Paul Eilers, Jan Gerretzen

References

Eilers, P.H.C. Eilers, P.H.C. (2004) "Parametric Time Warping", Analytical Chemistry, **76** (2), 404 – 411.

Boelens, H.F.M., Eilers, P.H.C., Hankemeier, T. (2005) "Sign constraints improve the detection of differences between complex spectral data sets: LC-IR as an example", Analytical Chemistry, **77**, 7998 – 8007.

Examples

data(gaschrom)
plot(gaschrom[1,], type = "1", ylim = c(0, 100))
lines(asysm(gaschrom[1,]), col = 2)

baseline.corr Baseline Correction using asymmetric least squares

Description

This function estimates a baseline using asymmetric least squares and subtracts it from the data.

Usage

baseline.corr(y, ...)

Arguments

У	signal(s) to correct. This can be a vector (containing one signal) or a matrix of signals(one signal per row)
	other arguments to the asysm function.

4

bestref

Value

ycorr

baseline corrected signal(s): a vector or a matrix of the same dimension as the input signal(s)

Author(s)

Paul Eilers, Jan Gerretzen

References

Eilers, P.H.C. (2004) "Parametric Time Warping", Analytical Chemistry, 76 (2), 404 - 411.

Examples

```
data(gaschrom)
plot(gaschrom[1,], type = "1", ylim = c(0, 100))
lines(baseline.corr(gaschrom[1,]), col = 2)
```

```
bestref
```

Identification of optimal reference

Description

This function calculates a similarity matrix and returns the sample number that is most similar to all other samples. This is possibly preferable as a reference sample since warping then may be kept to a minimum. Either RMS or WCC may be used as similarity functions.

Usage

Arguments

х	data matrix or array of signals, specified row-wise. In case of an array, the third dimension should differentiate between the different samples
optim.crit	either "WCC" or "RMS"
trwdth	the width of the triangle in the WCC criterion, given as a number of data points. Default: 20
wghts	Optional weights vector in the WCC criterion; will be calculated from the trian- gle width if necessary. Sometimes it is more efficient to pre-calculate it and give it as an argument
smooth.param	smoothing parameter for smoothing the signal when optim.crit equals "RMS". If no smoothing is required, set this to 0

Value

A list containing two elements:

best.ref	the index of the best reference(s)
crit.values	the qualities as measured by either RMS or WCC

Author(s)

Jan Gerretzen, Ron Wehrens

Examples

```
data(gaschrom)
bestref(gaschrom, optim.crit = "WCC", trwdth = 50)
bestref(gaschrom, optim.crit = "RMS")
bestref(gaschrom, optim.crit = "RMS", smooth.param = 1e5)
```

calc.multicoef	Calculation of warping coefficients when applying more than one
	warping function successively

Description

Applying two (or more) warping function after each other can be described with one warping function of a higher warping degree. This function provides the coefficients of this higher degree warping function.

Usage

```
calc.multicoef(coef1, coef2)
```

Arguments

coef1	vector containing the warping coefficients of the first applied warping function
coef2	vector containing the warping coefficients of the second applied warping func-
	tion

Details

This function uses Pascal's simplex to calculate the new warping coefficients.

When applying three warping functions successively (first a, then b and finally c - here a, b and c are vectors of warping coefficients), first calculate the new coefficients for b and c, and afterwards the coefficients for a with these new coefficients. So the coefficients for the total warping function can be calculated via calc.multicoef(a,calc.multicoef(b,c)).

6

calc.zerocoef

Value

a vector containing the corrected warping coefficients

Author(s)

Jan Gerretzen

References

Bloemberg, T.G., et al. (2010) "Improved parametric time warping for Proteomics", Chemometrics and Intelligent Laboratory Systems, 104 (1), 65 - 74.

See Also

calc.zerocoef

Examples

```
data(gaschrom)
ref <- gaschrom[1,]
samp <- gaschrom[16,]
coef1 <- c(100,1.1, 1e-5)
coef2 <- c(25, 0.95, 3.2e-5)
gaschrom.ptw <- ptw(ref, samp, init.coef = coef1, try = TRUE)
ref.w <- gaschrom.ptw$reference
samp.w <- gaschrom.ptw$warped.sample
samp.w[is.na(samp.w)] <- 0
gaschrom.ptw2 <- ptw(ref.w, samp.w, init.coef = coef2, try = TRUE)
plot(c(gaschrom.ptw2$warped.sample), type = "1")
corr.coef <- calc.multicoef(coef1, coef2)
gaschrom.ptw3 <- ptw(ref, samp, init.coef = corr.coef, try = TRUE)
lines(c(gaschrom.ptw3$warped.sample), col = 2, lty = 2)</pre>
```

calc.zerocoef Correction for warping coefficients when using zeropadding

Description

This function calculates the warping coefficients for the original range of the data, based on the warping of zero-filled data. Only needed when zeros are added in the beginning of the signal.

Usage

```
calc.zerocoef(coef, zeros)
```

Arguments

coef	vector of warping coefficients of a PTW-calculation on a set of signals with
	zeros added to the beginning of the signal
zeros	the number of zeros added

Value

a vector containing the corrected warping coefficients

Author(s)

Jan Gerretzen

References

Bloemberg, T.G., et al. (2010) "Improved parametric time warping for Proteomics", Chemometrics and Intelligent Laboratory Systems, **104** (1), 65 – 74.

See Also

padzeros calc.multicoef

Examples

```
data(gaschrom)
gaschrom.zf <- padzeros(gaschrom, 250)
ref <- gaschrom[1,]
samp <- gaschrom[16,]
ref.zf <- gaschrom.zf[1,]
samp.zf <- gaschrom.zf[16,]
gaschrom.ptw <- ptw(ref.zf, samp.zf)
layout(matrix(1:2,2,1, byrow=TRUE))
plot(gaschrom.ptw)
corr.coef <- calc.zerocoef(gaschrom.ptw$warp.coef, 250)
gaschrom.ptw2 <- ptw(ref, samp, init.coef = corr.coef, try = TRUE)
plot(gaschrom.ptw2)</pre>
```

coda

Chromatogram selection using the CODA algorithm

Description

The CODA algorithm calculates a so-called MCQ (Mass Chromatogram Quality) value for every row of the input. High MCQ values correspond with those chromatograms not containing spikes and/or a baseline.

Usage

```
coda(x, window = 5, smoothing = c("median", "mean"))
```

difsm

Arguments

х	data matrix containing chromatograms in the rows
window	width of the smoothing window
smoothing	type of smoothing: whether to use running means or running medians

Details

The MCQ value of a spectrum is the inner product between the standardized, smoothed chromatogram, and the length-scaled chromatogram. In literature, a cut-off of 0.85 has been reported to work well in selecting useful chromatograms, although this is strongly data-set dependent.

References

Windig, W., Phalp, J., Payna, A. (1996) "A noise and background reduction method for component detection in liquid chromatography/mass spectrometry", Analytical Chemistry, **68**, 3602 – 3606.

Examples

data(gaschrom)
coda(gaschrom)

difsm

Smoothing with a finite difference penalty

Description

This function smoothes signals with a finite difference penalty of order 2.

Usage

difsm(y, lambda)

Arguments

У	signal to be smoothed: a vector
lambda	smoothing parameter: larger values lead to smoothing

Value

smoothed signal: a vector

Author(s)

Paul Eilers, Jan Gerretzen

References

Eilers, P.H.C. (2004) "Parametric Time Warping", Analytical Chemistry, **76** (2), 404 – 411. Eilers, P.H.C. (2003) "A perfect smoother", Analytical Chemistry, **75**, 3631 – 3636.

Examples

```
data(gaschrom)
plot(gaschrom[1,], type = "1", ylim = c(0, 100))
lines(difsm(gaschrom[1,], lambda = 1e5), col = 2)
lines(difsm(gaschrom[1,], lambda = 1e6), col = 3)
lines(difsm(gaschrom[1,], lambda = 1e7), col = 4)
```

gaschrom

16 calibration GC traces

Description

The object gaschrom contains 16 calibration GC traces, measured at 5,000 time points. A peakpicked version is available as object gaschrom.st (see example).

Usage

data(gaschrom)

Source

Claire Boucon, Unilever

References

Eilers, P.H.C. (2004) "Parametric Time Warping", Analytical Chemistry, 76 (2), 404 - 411.

Examples

```
data(gaschrom)
```

```
## the gaschrom.st object has been generated in the following way:
## Not run:
pick.peaks <- function(x, span) {
   span.width <- span * 2 + 1
   loc.max <- span.width + 1 -
        apply(embed(x, span.width), 1, which.max)
   loc.max[loc.max == 1 | loc.max == span.width] <- NA
   pks <- loc.max + 0:(length(loc.max)-1)
   pks <- pks[!is.na(pks)]
   pks.tab <- table(pks)</pre>
```

lcms

Parts of 3 proteomic LC-MS samples

Description

The lcms data consists of a 100 x 2000 x 3 array lcms, a vector time of length 2000 and a vector mz of length 100. The LC-MS data in the array are a subset (samples 1, 2 and 5) of a larger set measured on a tryptic digest of E. coli proteins (see source section). Separate objects mz and rt give the values for the first two axis.

Peak picking leads to the object ldms.pks (see example section).

Usage

data(lcms)

Source

Nijmegen Proteomics Facility, Department of Laboratory Medicine, Radboud University Nijmegen Medical Centre. Data available (in different formats) at http://www.cac.science.ru.nl/research/data/ecoli/

References

Bloemberg, T.G., et al. (2010) "Improved parametric time warping for Proteomics", Chemometrics and Intelligent Laboratory Systems, **104** (1), 65 – 74.

Examples

```
## the lcms.pks object is generated in the following way:
## Not run:
data(lcms)
pick.peaks <- function(x, span) {</pre>
  span.width <- span * 2 + 1</pre>
  loc.max <- span.width + 1 -</pre>
      apply(embed(x, span.width), 1, which.max)
  loc.max[loc.max == 1 | loc.max == span.width] <- NA</pre>
  pks <- loc.max + 0:(length(loc.max)-1)</pre>
  pks <- pks[!is.na(pks)]</pre>
  pks.tab <- table(pks)</pre>
  pks.id <- as.numeric(names(pks.tab)[pks.tab > span])
  cbind(rt = pks.id, I = x[pks.id])
}
## bring all samples to the same scale, copied from ptw man page
lcms.scaled <- aperm(apply(lcms, c(1,3),</pre>
                             function(x) x/mean(x) ), c(2,1,3))
lcms.s.z <- aperm(apply(lcms.scaled, c(1,3),</pre>
                          function(x) padzeros(x, 250), c(2,1,3))
lcms.pks <- lapply(1:3,</pre>
                    function(ii) {
                      lapply(1:nrow(lcms.s.z[,,ii]),
                              function(jj)
                              cbind("mz" = jj,
                                    pick.peaks(lcms.s.z[jj,,ii], 5)))
                    })
## End(Not run)
```

mzchannel2pktab

Conversion between peak lists from hyphenated MS (LCMS, GCMS, ...) data and input for stptw.

Description

Function pktab2mzchannel allows to split a list of peaks into several sublists, for instance on the basis of m/z values. The result can be aligned with stptw. The peak list can be obtained from packages like XCMS. The reverse function, mzchannel2pktab, simply gathers all peak positions in one matrix.

Usage

```
pktab2mzchannel(pktab, Ivalue = "maxo", masses = NULL, nMasses = 0, massDigits = 2)
mzchannel2pktab(mzchannels)
```

12

padzeros

Arguments

pktab	a peak table as generated for example by XCMS. Necessary information: m/z value (column name "mz"), retention time (column name "rt") and intensity.
Ivalue	the name of the intensity value to be used. Default is "maxo", one of the columns generated by the XCMS package.
masses	a vector of specific masses
nMasses	an optional number limiting the number of mass channels. When both masses and nMasses are defined, the former takes preference
massDigits	number of significant mass digits - if no masses are supplied this number deter- mines the number of distinct categories in the output
mzchannels	a list of peak matrices, e.g. the output of pktab2mzchannel

Value

Function pktab2mzchannel returns a list of peak matrices; list elements have the name of the mz value that they represent. Function mzchannel2pktab returns one peak matrix where all masses are in a specific column.

Author(s)

Ron Wehrens

Examples

```
data(lcms)
## first couple of peaks in the first three channels
(smallset <- lapply(lcms.pks[[1]][1:3], head))
## all in one data matrix
allpeaks <- mzchannel2pktab(smallset)
## and back again
pktab2mzchannel(allpeaks, Ivalue = "I")</pre>
```

padzeros

```
Pad matrix with zeros
```

Description

Adds zeros to the left side of a matrix or vector, to its right side, or to both sides.

Usage

```
padzeros(data, nzeros, side="both")
```

Arguments

data	the original matrix or vector
nzeros	number of columns to add on one side
side	to which side to add the zeros - choose between "both", "left" or "right"

Details

When data is a numeric vector, it is converted to a matrix of a single row.

Value

A matrix with the same number of rows as the original matrix, and extra columns containing zeros on the specified side or sides

Author(s)

Jan Gerretzen

References

Bloemberg, T.G., et al. (2010) "Improved parametric time warping for Proteomics", Chemometrics and Intelligent Laboratory Systems, **104** (1), 65 – 74.

Examples

```
data(lcms)
lcms.z1 <- padzeros(lcms[75,,1], 250, side="left")
lcms.z2 <- padzeros(lcms[75,,1], 250, side="right")
lcms.z3 <- padzeros(lcms[75,,1], 250, side="both")
zeros <- rep(0, 250)
layout(matrix(1:4,2,2, byrow=TRUE))
plot(lcms[75,,1], type="1", main="Original signal")
plot(as.vector(lcms.z1), type="1", main="Padzeros left side")
points(1:250, zeros, col=2, lwd=0.08)
plot(as.vector(lcms.z2), type="1", main="Padzeros right side")
points(2001:2250, zeros, col=2, lwd=0.08)
plot(as.vector(lcms.z3), type="1", main="Padzeros both sides")
points(1:250, zeros, col=2, lwd=0.08)
```

points(2251:2500, zeros, col=2, lwd=0.08)

plot.ptw

Description

The function plots a ptw object. It shows either the original and warped signals, or the warping function.

Usage

Arguments

x	object of class 'ptw'
what	"signal" plots the reference, sample and warped sample signal(s); "function" plots the warping function as warped 'time' - 'time' for the forward warping mode and as 'time' - warped 'time' for the backward warping mode.
type	"simultaneous" plots all signals or warping functions in one plot; "individual" generates multiple plots
ask	logical, whether to ask before showing a new plot
	further arguments to the plotting functions

Author(s)

Jan Gerretzen, Ron Wehrens, Tom Bloemberg

Examples

predict.ptw

Description

Given a ptw object, predict either the signal at a certain warped time, or the warped time itself.

Usage

```
## S3 method for class 'ptw'
predict(object, newdata, what = c("response", "time"),
RTref = NULL, ...)
```

Arguments

object	An object of class "ptw"
newdata	Optional vector or matrix of new data points. If what equals "response", the new data should be a vector or matrix of intensities. If what equals "time", the new data is a vector of time points (a matrix of time points makes no sense).
what	Either "response", in which case the function returns the warped signal, or "time", and then the function returns the warped time axis. That is, the time point in the warped sample corresponding to the given time point in the original sample.
RTref	Optional vector of retention times in the reference.
	Further arguments, at the moment not used.

Value

The function returns a matrix (possibly containing only one row) of either warped time points or signals, warped according to the warping function defined in object. When warping signals individually, predict.ptw will check the dimension of newdata: if this is a vector or a matrix of one row, every single warping function will be applied to the one row. If the number of rows equals the number of warping functions, each row will be warped with its corresponding function. If the number of rows does not match the number of warping functions and is not equal to one, an error is given.

Author(s)

Ron Wehrens

References

Eilers, P.H.C. "Parametric Time Warping." Anal. Chem., 2004, 76, 404-411 Bloemberg, T.G. et al. "Improved parametric time warping for proteomics." Chemom. Intell. Lab. Syst., 2010, 104, pp. 65-74

predict.ptw

See Also

ptw

Examples

```
## educational example, contributed by zeehio (Sergio Oller)
x1 <- c(rep(0, 5), 1,1,1, 20, 40, 20, 1, 1, 1, rep(0, 5))
x2 <- c(rep(0, 6), 1,1,1, 20, 40, 20, 1, 1, 1, rep(0, 4))
time <- 1:length(x1)
## get time-warped object. Default: 'forward' warping, also works
## with backward warping
w1b <- ptw(ref = x1, samp = x2)
## predict intensities of object x2 after warping at the times used in x1
x2wb <- predict(w1b, newdata = x2, what = "response")
## predict times where the original elements of x2 will end up
t2wb <- as.numeric(predict(w1b, newdata = time, what = "time"))
graphics.off()
par(mfrow = c(2,1))
plot(x1, type = "h", col = 2, lwd = 2, main = "Orig data")
points(x2, type = "h", col = 4)
plot(x1, type = "h", col = 2, lwd = 2, main = "Backward warping")
points(c(x2wb), type = "h", col = 4) # what = "response"
points(t2wb, x2, col = 4)
                                    # what = "time"
## more relevant example
data(gaschrom)
## Global warping: all samples warped with the same function
ref <- gaschrom[1,]</pre>
samp <- gaschrom[14:16,]</pre>
gp <- ptw(ref, samp, init.coef = c(0, 1), warp.type = "global")</pre>
matplot(t(samp), type = "1", xlim = c(2200, 2400), lty = 1, col = 1:3)
lines(ref, type = "1", col = "gray", lwd = 2)
## plot predicted warped signal directly
matlines(t(predict(gp)), lty = 2, col = 1:3)
## plot original signal at warped time axis
matlines(t(predict(gp, newdata = 2001:2600, what = "time")),
         t(samp[,2001:2600]), col = 1:3, lwd = 3, lty = 2) ## OK
## result: good alignment with ref, differences between three profiles persist
## Individual warping: all samples warped individually
gp <- ptw(ref, samp, init.coef = c(0, 1), warp.type = "indiv")</pre>
predict(gp, what = "time", newdata = 2001:2600)
matplot(t(samp), type = "1", xlim = c(2200, 2400), lty = 1, col = 1:3)
lines(ref, type = "1", col = "gray", lwd = 2)
matlines(t(predict(gp, what = "time")),
         t(samp), col = 1:3, lty = 2)
## result: each individual profile is aligned to the ref
## How would samples 11:13 be warped using the coefficients from samples
## 14:16 (silly but just to make the point)?
```

samp.pred <- predict(gp, what = "response", newdata = gaschrom[11:13,])</pre>

ptw

Parametric Time Warping

Description

The main functions of the ptw package, performing parametric time warping of one or more samples. Features in the samples are optimally aligned with features in the reference(s). One may align a single sample to a single reference, several samples to a single reference, and several samples to several references. In the latter case, the number of references and samples should be equal. One may require that all samples are warped with the same warping function, or one may allow individual warpings for all samples.

Usage

```
ptw(ref, samp, selected.traces,
    init.coef = c(0, 1, 0), try = FALSE,
    warp.type = c("individual", "global"),
    optim.crit = c("WCC", "RMS"),
    mode = c("forward", "backward"),
    smooth.param = ifelse(try, 0, 1e05),
    trwdth = 20, trwdth.res = trwdth,
    verbose = FALSE, ...)
stptw(ref, samp, init.coef = c(0, 1, 0),
    trwdth = 20, trwdth.res = trwdth,
    nGlobal = ifelse(length(init.coef) > 3, 5, 0),
    ...)
## S3 method for class 'ptw'
summary(object, ...)
## S3 method for class 'ptw'
print(x, ...)
```

Arguments

ref	reference. Either a vector (containing one reference signal) or a matrix (one reference per row). If more than one reference is specified, the number of reference signals must equal the number of sample signals.	
samp	sample. A vector (containing one sample signal) or a matrix (one sample per row).	
selected.traces		
	optional vector containing the row numbers to use from ref (if more than one reference signal is specified) and samp.	
init.coef	starting coefficients. The first number is the zeroth-order coefficient (i.e., a con- stant shift); further numbers indicate linear, quadratic, stretches. The default is to start from the identity warping using a quadratic function $(c(0,1,0))$	

try

warp.type

optim.crit

mode

if try = TRUE, ptw does not optimize the warping but returns a ptw object con- taining the warping for init.coef. Default: FALSE
default is to treat samples and references as single entities and align them indi- vidually and independently. Using the argument warp.type = "global" leads to one alignment function; the samples are warped simultaneously to the refer- ence(s). Also see details
either "WCC" or "RMS". In both cases, the optimal value of the alignment leads to a value of 0. For "WCC", this means that 1 -WCC is optimized rather than WCC (where the optimal value equals 1)
either "forward" or "backward": the latter was the original implementation, ba- sically for a point i in the original signal predicting the point j in the signal that would be in position i in the warped signal. The interpretation of the coefficients is counterintuitive. Therefore the default is "forward", simply predicting the lo- cation (time) in the warped signal of a particular point. Apart from possible numerical optimisation issues, both warpings should give the same net result.
smoothing parameter for smoothing the reference and sample when optim.crit

smooth.param	smoothing parameter for smoothing the reference and sample when optim.crit equals "RMS". If no smoothing is required, set this to 0. The default is to use smoothing in the optimization mode, and no smoothing otherwise
trwdth	the width of the triangle in the WCC criterion during the optimization, given as a number of data points. Default: 20
trwdth.res	the width of the triangle in the WCC calculation in the calculation of the quality of the final result. Default: equal to $trwdth$
verbose	logical, default is FALSE. Whether to give output during the optimisation, which may be useful for large data sets
	further arguments to optim
nGlobal	for stick optimisations with more than three coefficients, a more powerful opti- misation method is used (nloptr: in such cases nGlobal indicates the number of times this method is employed.

x, object an object of class "ptw"

Details

Function ptw and friends is meant for profile data, where intensities have been recorded at regular time points; function stptw is meant for lists of peaks, for instance obtained after peak-picking the profile data. The latter option is less flexible (Euclidean distance and backward warping have not been implemented, and only global alignment is possible) but is much faster, especially for larger data sets.

In the optimization mode (try = FALSE), the function optimizes the warping coefficients using the chosen criterion (either "WCC" or "RMS"). For "RMS", the data are smoothed before the optimization, but the quality of the final warping is measured on the unsmoothed data. For "WCC", the warping is performed using trwdth as the triangle width, but the quality of the final solution is measured using trwdth.res.

If try = TRUE is used as an argument, the function does not start an optimization, but just calculates the warping for the given warp function (init.coef); if smooth.param is larger than zero for the RMS criterion, the RMS of the smoothed patterns is calculated. The WCC criterion uses trwidth.res as the triangle width in this case.

Five situations can be distinguished:

- 1. One sample and one reference: this obviously leads to one warping function regardless of the setting of warp.type.
- 2. Several samples, all warped to the same single reference, each with its own warping function: this is the default behaviour (warp.type = "individual")
- 3. Several samples, warped to an equal number of references (pair-wise), with their own warping functions: this is the default behaviour (warp.type = "individual")
- 4. Several samples, warped to one reference, with one warping function (warp.type = "global")
- Several samples, warped to an equal number of references (pair-wise), with one warping function (warp.type = "global")

Value

A list of class "ptw" containing:

reference	the reference(s) used as input
sample	the sample(s) used as input
warped.sample	the warped sample
warp.coef	the warping coefficients
warp.fun	the warped indices (not for stptw)
crit.value	the value of the chosen criterion, either "WCC" or "RMS"
optim.crit	the chosen criterion, either "WCC" or "RMS"
warp.type	the chosen type of warping, either "individual" or "global"

Author(s)

Jan Gerretzen, Ron Wehrens

References

Eilers, P.H.C. (2004) "Parametric Time Warping", Analytical Chemistry, 76 (2), 404 - 411.

Bloemberg, T.G., et al. (2010) "Improved parametric time warping for Proteomics", Chemometrics and Intelligent Laboratory Systems, 104 (1), 65 - 74.

See Also

WCC, RMS, select.traces, gaschrom, lcms

ptw

Examples

```
data(gaschrom)
ref <- gaschrom[1,]</pre>
samp <- gaschrom[16,]</pre>
gaschrom.ptw <- ptw(ref, samp)</pre>
summary(gaschrom.ptw)
## same with sticks (peak lists)
refst <- gaschrom.st[1]</pre>
sampst <- gaschrom.st[16]</pre>
gaschrom.st.ptw <- stptw(refst, sampst, trwdth = 100)</pre>
summary(gaschrom.st.ptw)
## Not run:
## comparison between backward and forward warping
gaschrom.ptw <- ptw(ref, samp, init.coef = c(0, 1, 0, 0), mode = "backward")</pre>
summary(gaschrom.ptw)
gaschrom.ptw \langle - ptw(ref, samp, init.coef = c(-10, 1, 0, 0), mode = "forward")
summary(gaschrom.ptw)
## many samples warped on one reference
ref <- gaschrom[1,]</pre>
samp <- gaschrom[2:16,]</pre>
gaschrom.ptw <-
 ptw(ref, samp, warp.type = "individual", verbose = TRUE,
      optim.crit = "WCC", trwdth = 100, init.coef = c(0, 1, 0))
summary(gaschrom.ptw)
## "individual" warping not implemented for sticks; do separate warpings
## instead
refst <- gaschrom.st[1]</pre>
sampst <- gaschrom.st[2:16]</pre>
gaschrom.st.ptw.list <-</pre>
 lapply(sampst,
         function(smpl)
           stptw(refst, list(smpl), trwdth = 100, init.coef = c(0, 1, 0)))
t(sapply(gaschrom.st.ptw.list, "[[", "warp.coef"))
t(sapply(gaschrom.st.ptw.list, "[[", "crit.value"))
## several samples on several references individually
ref <- gaschrom[1:8,]</pre>
samp <- gaschrom[9:16,]</pre>
gaschrom.ptw <-
 ptw(ref, samp, warp.type = "individual",
      optim.crit = "WCC", trwdth = 100, init.coef = c(0, 1, 0))
summary(gaschrom.ptw)
## stick version
```

gaschrom.st.ptw.list <mapply(function(x, y)</pre>

```
stptw(list(x), list(y), trwdth = 100, init.coef = c(0, 1, 0)),
        gaschrom.st[1:8], gaschrom.st[9:16],
        SIMPLIFY = FALSE)
t(sapply(gaschrom.st.ptw.list, coef))
## several samples on several references: one, global warping
gaschrom.ptw <- ptw(ref, samp, warp.type = "global",</pre>
                   optim.crit = "WCC", init.coef = c(0, 1, 0))
summary(gaschrom.ptw)
refst <- gaschrom.st[1:8]</pre>
sampst <- gaschrom.st[9:16]</pre>
gaschrom.st.ptw <- stptw(refst, sampst, trwdth=100, init.coef = c(0, 1, 0))</pre>
summary(gaschrom.st.ptw)
** ******************
## Example of a three-way data set#
** *****************
## first bring all samples to the same scale
data(lcms)
lcms.scaled <- aperm(apply(lcms, c(1,3),</pre>
                         function(x) x/mean(x) ), c(2,1,3))
## add zeros to the start and end of the chromatograms
lcms.s.z <- aperm(apply(lcms.scaled, c(1,3),</pre>
                       function(x) padzeros(x, 250) ), c(2,1,3))
## define a global 2nd degree warping
warp1 <- ptw(lcms.s.z[,,2], lcms.s.z[,,3], warp.type="global")</pre>
warp.samp <- warp1$warped.sample</pre>
warp.samp[is.na(warp.samp)] <- 0</pre>
## refine by adding 5th degree warpings for individual chromatograms
warp2 <- ptw(lcms.s.z[,,2], warp.samp, init.coef=c(0,1,0,0,0,0))</pre>
warp.samp2 <- warp2$warped.sample</pre>
warp.samp2[is.na(warp.samp2)] <- 0</pre>
## compare TICs
layout(matrix(1:2,2,1, byrow=TRUE))
plot(colSums(lcms.s.z[,,2]), type="1", ylab = "",
    main = "TIC: original data")
lines(colSums(lcms.s.z[,,3]), col=2, lty=2)
plot(colSums(lcms.s.z[,,2]), type="l", ylab = "",
    main = "TIC: warped data")
lines(colSums(warp.samp2), lty=2, col=2)
## stick version of this warping - note that the peaks have been picked
## from the scaled profiles. Note that here we need to take list
## elements: every sample is a list of mz channels.
```

warp1.st <- stptw(lcms.pks[[2]], lcms.pks[[3]], trwdth = 100) summary(warp1.st)

ptwgrid

End(Not run)

ptwgrid

Calculate RMS or WCC values on a grid

Description

Calculates values of the chosen optimization criterion (RMS or WCC) on a grid defined by the coefficients of the warping function.

Usage

```
ptwgrid(ref, samp, selected.traces,
coef.mins, coef.maxs, coef.lengths,
optim.crit = c("WCC", "RMS"),
smooth.param = 1e05,
trwdth = 20)
```

Arguments

ref	reference. Either a vector (containing one reference signal) or a matrix (one ref- erence per row). If more than one reference is specified, the number of reference signals must equal the number of sample signals
samp	sample. A vector (containing one sample signal) or a matrix (one sample per row). If more than one sample is specified, the number of sample signals must equal the number of reference signals
selected.traces	
	optional vector containing the row numbers to use from ref (if more than one reference signal is specified) and samp
coef.mins	a vector containing the respective minimal values of coefficients for the grid. The first number is the minimal zeroth-order coefficient (i.e., a constant shift); further numbers indicate the minimal linear, quadratic, stretches
coef.maxs	a vector containing the maximal values of coefficients for the grid
coef.lengths	a vector of the same length as coef.maxs and coef.mins containing the number of steps in which to vary the respective coefficients between their minimum and maximum value
optim.crit	either "WCC" or "RMS". In both cases, the optimal value of the alignment leads to a value of 0. For "WCC", this means that 1 -WCC is optimized rather than WCC (where the optimal value equals 1)
smooth.param	smoothing parameter for smoothing the reference and sample when optim.crit equals "RMS". If no smoothing is required, set this to 0.
trwdth	the width of the triangle in the WCC criterion during the optimization, given as a number of data points. Default: 20

Details

In contrast to ptw, only the three situations leading to one warping function are distinguished here:

- 1. One sample and one reference;
- 2. Several samples, warped to an equal number of references (pair-wise);
- 3. Several samples, warped to a single reference.

Which situation is applicable is determined from the dimensions of ref and samp.

Value

An array of dimension length(coef.mins) and maximal indices per dimension specified by coef.lengths

Author(s)

Tom Bloemberg, Jan Gerretzen, Ron Wehrens

See Also

ptw

Examples

End(Not run)

RMS

Quality criteria for comparing patterns with shifts

Description

Functions to compare patterns with shifted features. These functions compare warped sample patterns to one or more reference patterns. RMS returns the usual root-mean-squared difference measure; WCC returns 1-wcc, where wcc indicates the weighted cross-correlation. Perfect alignment leads to a value of 0 for both criteria.

Internal function, not meant to be called directly by the user. In particular, note that the identity warping may lead to slightly different estimates than a direct comparison of the reference and sample signals - a warping even slightly outside the original range of 1 : ncol(ref) leads to NA values.

RMS

Usage

```
RMS(warp.coef, ref, samp, B, mode)
WCC(warp.coef, ref, samp, B, trwdth = 20, wghts, mode, ref.acors = NULL)
```

Arguments

warp.coef	a vector of warping coefficients
ref	reference signal; a matrix with one or more rows. If the number of rows is greater than one, it should equal the number of rows in samp
samp	sample signal; a matrix with one or more rows
В	basis for warping function
mode	either "forward" (new implementation, also used for warping peak lists) or "back-ward" (classical implementation).
trwdth	triangle width for the WCC function, expressed in the number of data points
wghts	optional weights vector, will be calculated from triangle width if necessary. Sometimes it is more efficient to pre-calculate it and give it as an argument
ref.acors	autocorrelation of the reference. Since the reference is often unchanged over multiple evaluations (e.g., during an optimization), it is useful to pre-calculate this number

Details

All patterns in samp are warped using the same warping function, and then compared to ref, either pair-wise (when ref and samp are of the same size), or with the one pattern in ref.

Value

One number - either the WCC or RMS value

Author(s)

Jan Gerretzen, Tom Bloemberg, Ron Wehrens

References

Eilers, P.H.C. (2004) "Parametric Time Warping", Analytical Chemistry, 76 (2), 404 - 411.

de Gelder, R., Wehrens, R. and Hageman, J.A. (2001) "A generalized expression for the similarity of spectra: Application to powder diffraction pattern classification", Journal of Computational Chemistry, **22**, 273 – 289.

See Also

wcc

```
select.traces
```

Description

For alignment purposes, it may be useful to select traces which show clear features, and to throw away traces that contain mainly noise. This function implements three ways to achieve this: CODA, a criterion similar to varimax, and a criterion based on the highest intensity.

Usage

Arguments

X	a data matrix or an array. The first dimension signifies the traces from which a selection is to be made. If X is a matrix, the first usually corresponds to samples and the second dimension is the spectral dimension. If X is an array, the data are assumed to come from a hyphenated experiment, with the first dimension the chromatographic dimension, the second the spectral dimension and the third dimension corresponding to samples	
criterion	either Windig's CODA algorithm, a criterion calculating the variances of the length-scaled spectra, or a criterion giving the height of the highest peak	
window, smoothing		
	arguments to the coda function.	

Details

The CODA criterion in essence selects traces with no baseline and no spikes, but still containing significant intensities. The variance criterion aims at something similar: it calculates the variance (or standard deviation) of every trace after length scaling - traces with a high value show few highly structured features, whereas traces with a low value show noise or a significant baseline. The intensity criterion simply returns the intensity of the highest peak. The latter two criteria are simpler than CODA but implicitly assume that the traces have been preprocessed (i.c., spikes have been removed).

Value

The function returns a list with components

crit.vala vector containing the values of the criterion for all traces. If X is an array, the
function is recursively applied to all samples (elements of the third dimension)
- the results are multiplied to obtain one criterion value per trace

trace.nrs the order of the traces (from large to small)

warp.time

Author(s)

Ron Wehrens

See Also

coda

Examples

```
data(lcms)
ntrace <- dim(lcms)[1]
lcms.selection <- select.traces(lcms[,,1:2], criterion = "var")
good <- lcms.selection$trace.nrs[1]
bad <- lcms.selection$trace.nrs[ntrace]
par(mfrow = c(1,2))
matplot(lcms[good,,1:2], type = 'l', lty = 1)
matplot(lcms[bad,,1:2], type = 'l', lty = 1)</pre>
```

```
warp.time
```

Transform time according to a given warping function

Description

Warp time points according to a warping function.

Usage

warp.time(tp, coef)

Arguments

tp	A vector of time points, not necessarily equidistant.
coef	The coefficients of the parametric time warping function: the first coefficient is
	the constant shift, the second the linear stretch etcetera.

Value

The function returns a vector of warped time points.

Author(s)

Ron Wehrens

References

Wehrens, R. et al. (2015) "Fast parametric warping of peak lists", Bioinformatics. DOI: 10.1093/bioinformatics/btv299.

Examples

```
time <- 1:100
## simple shift and compression
warp.time(time, c(-10, .99))</pre>
```

WCC

Weighted auto- and cross-correlation measures

Description

Functions to calculate weighted auto- and cross-correlation measures. The wcc is a suitable measure for the similarity of two patterns when features may be shifted. Identical patterns lead to a wcc value of 1.

Functions wcc and wac are meant for profile data (intensities measured at equidistant time points), whereas wcc.st and wac.st are meant for peak lists. In general, wcc values calculated for profiles will be higher since they will also include the large similarity in the empty spaces, i.e., parts of the profiles where no peaks are present (and that will appear to be perfectly aligned), whereas the peak-based version concentrates only on the peaks.

Usage

wcc(pattern1, pattern2, trwdth, wghts = NULL, acors1 = NULL, acors2 = NULL)
wac(pattern1, trwdth, wghts = NULL)

Arguments

pattern1, patte	rn2
	input patterns, typically spectra. Vectors
trwdth	triangle width, given in the number of data points for the profile functions, and in the actual retention times for the stick-based warpings.
wghts	optional weights vector, will be calculated from triangle width if necessary. Sometimes it is more efficient to pre-calculate it and give it as an argument
acors1, acors2	autocorrelations of the input patterns. If not provided, they are calculated

Details

Functions wcc and wac are defined such that the triangle width stands for the number of points on one side of **and including** the current point. Thus, a trwdth of 0 signifies a non-existent triangle and results in an error; a trwdth equal to 1 only includes the current point with weight 1 and no neighbouring points. For the stick-based equivalents, the units of the time axis are used for the triangle width.

Value

One number, the weighted autocorrelation or crosscorrelation

28

whit1

Author(s)

Ron Wehrens

References

de Gelder, R., Wehrens, R. and Hageman, J.A. (2001) "A generalized expression for the similarity of spectra: Application to powder diffraction pattern classification", Journal of Computational Chemistry, **22**, 273 – 289.

Examples

```
data(gaschrom)
wcc(gaschrom[1,], gaschrom[2,], 20)
```

```
wcc.st(gaschrom.st[[1]], gaschrom.st[[2]], 20)
```

whit1	Weighted	Whittaker	smoothing	with a	a first	order	finite	difference
	penalty							

Description

This function smoothes signals with a finite difference penalty of order 1.

Usage

whit1(y, lambda, w)

Arguments

У	signal to be smoothed: a vector
lambda	smoothing parameter: larger values lead to more smoothing
W	weights: a vector of same length as y. Default weights are equal to one

Value

smoothed signal: a vector

Author(s)

Paul Eilers, Jan Gerretzen

References

Eilers, P.H.C. (2004) "Parametric Time Warping", Analytical Chemistry, **76** (2), 404 – 411. Eilers, P.H.C. (2003) "A perfect smoother", Analytical Chemistry, **75**, 3631 – 3636.

Examples

```
data(gaschrom)
plot(gaschrom[1,], type = "l", ylim = c(0, 100))
lines(whit1(gaschrom[1,], lambda = 1e1), col = 2)
lines(whit1(gaschrom[1,], lambda = 1e2), col = 3)
lines(whit1(gaschrom[1,], lambda = 1e3), col = 4)
```

		٠		\sim
141	n	п.	+	,
VV		т.	L	~
		+	L	~

Weighted Whittaker smoothing with a second order finite difference penalty

Description

This function smoothes signals with a finite difference penalty of order 2.

Usage

whit2(y, lambda, w)

Arguments

У	signal to be smoothed: a vector
lambda	smoothing parameter: larger values lead to more smoothing
W	weights: a vector of same length as y. Default weights are equal to one

Value

smoothed signal: a vector

Author(s)

Paul Eilers, Jan Gerretzen

References

Eilers, P.H.C. (2004) "Parametric Time Warping", Analytical Chemistry, **76** (2), 404 – 411. Eilers, P.H.C. (2003) "A perfect smoother", Analytical Chemistry, **75**, 3631 – 3636.

Examples

```
data(gaschrom)
plot(gaschrom[1,], type = "1", ylim = c(0, 100))
lines(whit2(gaschrom[1,], lambda = 1e5), col = 2)
lines(whit2(gaschrom[1,], lambda = 1e6), col = 3)
lines(whit2(gaschrom[1,], lambda = 1e7), col = 4)
```

30

Index

```
* datasets
    gaschrom, 10
    lcms, 11
* manip
    asysm, 3
    baseline.corr,4
    bestref, 5
    calc.multicoef, 6
    calc.zerocoef,7
    coda, 8
    difsm, 9
    mzchannel2pktab, 12
    padzeros, 13
    plot.ptw, 15
    predict.ptw, 16
    ptw, 18
    ptwgrid, 23
    RMS, 24
    select.traces, 26
    warp.time, 27
    wcc, 28
    whit1,29
    whit2, 30
* package
    ptw-package, 2
asysm,3
baseline.corr,4
bestref, 5
calc.multicoef, 6, 8
calc.zerocoef, 7,7
coda, 8, 27
coef.ptw(ptw), 18
difsm, 9
gaschrom, 10, 20
lcms, 11, 20
```

mz (lcms), 11 mzchannel2pktab, 12 padzeros, 8, 13 pktab2mzchannel(mzchannel2pktab), 12 plot.ptw, 15 predict.ptw, 16 print.ptw(ptw), 18 print.stptw(ptw), 18 ptw, 17, 18, 24 ptw-package, 2 ptwgrid, 23 RMS, 20, 24 select.traces, 20, 26 stptw (ptw), 18 summary.ptw(ptw), 18 summary.stptw(ptw), 18 time (lcms), 11 wac (wcc), 28 warp.time, 27 WCC, 20 WCC (RMS), 24 wcc, 25, 28 whit1,29 whit2, 30